


Accessories

Accessories

ColorFlex® EZ



	Port Inserts
04-6623-10	32 mm (1.25 in)
02-6624-05	32 mm (1.25 in) with glass
D02-1010-618	32 mm (1.25 in) with 420nm UV Filter Assembly
04-6623-09	25 mm (1 in)
02-6624-01	25 mm (1 in) with glass
04-6623-08	19 mm (0.75 in)
04-6623-05	13 mm (0.50 in)



04-6622-00

Port Insert for Sample Cup

64 mm (2.5 in) The use of the sample cup enables the measurement of powders, granules, pellets, pastes and liquids. The port insert is designed to hold the sample cup in proper position at the measuring port.



04-7209-00

Glass Sample Cup: 64 mm (2.5 in)

Optically clear glass cup for sample presentation of liquids, powders, granules, and pellets.



02-4579-00

Ring and Disk Set

For use in the glass sample cup, used for transparent and translucent liquids on semi-solids where path length must be fixed. Does not include sample cup.



04-4000-00

Sample Cup Opaque Cover

Provides a light trap to exclude the interference of external light on the sample cup.



D02-1010-327

Footswitch Assembly

The footswitch is used to initiate sample measurement without the use of hands.

D02-1010-459

Sample Clamp / Port Forward Stand

Securely holds the instrument in a port forward orientation. The sample clamp permits to be held flat against the sample port.



L02-1014-582

Hard-Sided Carrying Case, CFEZ

Hard-sided Carrying Case for Standard CFEZ

L02-1014-607

Hard-Sided Carrying Case, CFEZ Citrus

Hard-sided Carrying Case for Citrus CFEZ

L02-1014-608

Hard-Sided Carrying Case, CFEZ Tomato

Hard-sided Carrying Case for Standard CFEZ



A13-1014-294

USB Flexible Keyboard

88 character keyboard allows user to enter IDs directly into CFEZ.



A13-1014-254

USB Barcode Reader

Barcode Scanner scans product IDs directly into CFEZ. Automatically detected at USB port.



A13-1014-259

USB Printer

USB printer allows hard copy measurements to be printed. requires a standard USB cable, sold separately. Uses Thermal Printer paper, sold separately.

A21-1013-859

USB Cable

USB cable for use with CFEZ and USB Printer and CFEZ.

L02-1014-440

Thermal Printer Paper

Thermal printer paper for use with A13-1014-259 USB printer. 5 rolls per package.

For more information, please contact HunterLab at 703-471-6870, sales@hunterlab.com or visit www.hunterlab.com